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Dental Bleaching Gel for Vital Teeth

Professional Use Only

Read all the information in this instruction manual before using the product. Save it for consultation, at least until the product has been totally used or there is no more interaction of the product with your last patient.

Description

The formulation of the dental bleaching gel Whiteness is based on products already known in the world market and equally combining the high viscosity of carbopol gel with the bleaching effect of carbamide peroxide to obtain effective dental bleaching. The pH of the gel is neutral so that it does not cause decalcification of the enamel and dentin. Whiteness Perfect is the bleaching gel which has, besides glycol and a high water content, potassium nitrate and sodium fluoride as desensitizers. The contents of the kit, if administered correctly, are sufficient for a complete and durable bleaching of both arches. The Mini kit contains enough gel for bleaching an arch or executing a complete retouch. The gel is applied in an individual tray. The gel concentrations 10% and 16% can be used during the night while the patient sleeps or during the day, depending on the dental surgeon's orientation. It is recommended for Whiteness Perfect 10% e 16% to be in contact with the teeth for a period of 3 to 4 hours daily (the gel's action period). Whiteness Perfect 22% can be used during 1 hour daily. Alternatively, longer or shorter times can be indicated according to the case.

In normal cases and when used daily, the treatment lasts from 10 to 14 days. The frequency or the time of using the gel can, however, be reduced to minimize sensitivity effects. Generally, the bleaching effects become visible after 3 to 5 days of use.

The Choice of Concentration

Whiteness Perfect is available in the following concentrations: 10%, 16% and 22%. In the cases of normal (natural) pigmentation by age or tobacco, coffee or also when there is a greater dental sensitivity the use of gel at 10% is indicated. In cases of greater pigmentation, more difficulty in removal or if visible results are desired in a shorter time, is indicated the use of gel with higher concentration at 16% or 22%.

In the cases of staining by tetracycline some improvement can be expected, but only in milder cases more satisfactory results can be obtained. Longer treatment times are generally required.

Forms of Presentation

Whiteness Perfect 10% and 16% Kit: 5 syringes of bleaching gel containing 3g each, 5 gel application tips, 2 vinyl plates 1mm thick for making the trays, 1 case for keeping the trays, and instructions for the professional and patient.

Kit Whiteness Perfect 22%: 4 syringes of bleaching gel containing 3g each, 4 gel application tips, 2 vinyl plates 1mm thick for making the trays, 1 case for keeping the trays, and instructions for the professional and patient.

Whiteness Perfect 10 and 16% Mini kit: 3 syringes of bleaching gel containing 3g each, 3 gel application tips, and instructions for the professional and patient.

Bulk: 50 syringes each with 3g of gel, packed individually.

Basic Composition

Carbamide Peroxide, Neutralized Carbopol, Potassium Nitrate, Sodium Fluoride, Humectant (Glycol), Deionized Water.

Indication

Whiteness Perfect is indicated for bleaching vital teeth using the technique of bleaching at home supervised by a dental surgeon. The product bleaches teeth darkened by age and naturally dark or yellow, and teeth stained by external pigments, such as coffee, tobacco, tea, etc.

Precautions and Contraindications

The use of the product must be indicated and always accompanied by a dental surgeon.

The oral health must be evaluated prior to treatment by means of clinical and radiographic examination.

Cracks in the enamel, defective restorations, areas of carie exposed root surfaces (may experience sensitivity) and other factors which can endanger the bleaching treatment must be restored before starting the treatment. The product has contraindications for patients who cannot solve the aforementioned problems.

The product is not indicated for bleaching teeth with amelogenesis and dentinogenesis imperfecta, severe fluorosis, intense stains by tetracycline and other anomalies of enamel and dentin which put at risk the vitality of the tooth.

As a safety precaution, the product should NOT be used during pregnancy or nursing period and under 15 years old.

For diabetic patients, it is recommended that the disease is properly controlled before and during the use of the bleaching gel in order to avoid irritant effects of the product. Alternatively, consider the use of an in-office bleaching gel which has no contact to gingival tissue.

The product is not recommended for people who have recently had a periodontal surgery or gingivitis or in case of the patient has a known allergy or chemical sensitivity to peroxides, carbopol, etc.

Restorations should be adequately sealed and all exposed sensitive dentin should be covered. If a history of sensitivity exists, treat with appropriate restoration or temporarily with dentin sealant.

Tell the patient about the level of results which can be obtained with the bleaching. Inform patients that existing crowns, tooth colored fillings, etc., do not whiten and will have to be changed at the end of the treatment (usually 7 days after the bleaching), teeth which are not vital can be bleached less intensely and that a special bleaching (Whiteness Super Endo) of this element may be required.

Try to obtain a photograph before and after the bleaching in order to see the results better. It is also a good idea to bleach one arch at a time. This will allow one arch to be a reference for the other one.

Tetracycline stains will require more treatment time than other types of disfiguring stains. Teeth with grey colouring, even though it is not caused by using tetracycline, may also not respond well to bleaching.

Occasionally a slight purple discoloration may appear in the tray following bleaching adjacent to amalgam fillings.

Teach the patient how to place the gel in the tray. We recommend one drop of gel in each tooth space of the tray. The excess of gel can be avoided and removed by using a toothbrush or his finger. Prolonged contact with the gingiva can result in localized irritation.

Give the patient the leaflet with the instructions and emphasize the importance of reading it.

Tell the patient about possible sensitivity to temperature variances, mainly to cold.

Exposed root surfaces may experience sensitivity or another adverse factor, protect the area or cut the tray to the cement union in order to avoid prolonged contact of the gel with the sensitive area. In these cases one cannot expect the uncovered area to bleach with the same intensity. There are cases where the sensitivity is not easy to treat. If none of the treatment forms described have an effect, using the gel can be prescribed for alternating days, every 3 days or for 1 to 2 hours a day prolonging the bleaching. In these cases it is advisable to use the gel with the lower concentration (10%). Generally the sensitivity disappears in at most 3 days. If the patient has too much discomfort or an allergic reaction to the product, the treatment must be suspended.

Guide the patient about not ingesting acid drinks or food as they can increase the sensitivity during the bleaching. Strongly coloured drinks and food should also be avoided so as not to affect the bleaching process.

Try to follow up the bleaching by asking the patient to return every 5 to 7 days, for example, in order to evaluate the treatment, mainly at the start of the bleaching.

Gingival irritation can also occur when using products for rinse the mouth or even toothpaste with a strong taste before using the tray with the gel. If you wish to use these products, it is advisable to leave an interval of at least 30 minutes (or sufficient time for the sensitivity to disappear) between using the mouthwash and or toothpaste and the gel. Aggressive brushing on the gingiva can also make it sensitive and result in irritation caused by the gel. Gingival irritation can also be caused by a custom tray which was not adequately trimmed.

Side Effects

Occasionally sensitivity may require treatment. If the sensitivity becomes uncomfortable, treat the teeth with the gel Desensibilize KF 0,2% (moderate to low sensitivity) or Desensibilize KF 2% (moderate to high sensitivity).

Areas of hypocalcification can exist which may not be visible to the naked eye. Hypocalcified areas will whiten faster, thereby becoming more obvious during bleaching. In some cases these stains disappear owing to the dehydration of the enamel after concluding the bleaching. Treating these areas with fluoride can also help to make them disappear.

Some patients can feel slight irritation in the gingiva, tongue or lips, generally caused by excessive use of gel in the tray.

Instructions for Use

See the instructions in the Patient's Instructions Manual.

Preservation and Storage

Store the product at 5°C and 25°C / 41°F and 77°F. Do not freeze the product. Keep it out of direct sunlight.

Warning

Do not use the product after expiry date. In order to discard the product, follow the legislation of your country. Maintain it out of reach of children.

Irritant to eyes and skin for concentrations between 5% to 8% of Hydrogen Peroxide (16% to 25.6% of Carbamide Peroxide).

This material was only manufactured for dental use and must be handled according to the instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Besides this, the user is obliged to prove, before use and under his or her responsibility that this material is compatible with the use desired, mainly when this use is not indicated in the instructions for use. Data descriptions do not constitute any type of guarantee and so are not in way linked.

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